

Product Specification



Product Codes	NPF 5001-5006
Product type	Examination and protective
UDI DI Code	Printed on the box and carton
Chlorinated or Not	Chlorinated
Color	White
Medical grade : (EU) 2017/745	Class I
Protection type: (EU) 2016/425	Category III
Food Contact:	Pass done by SGS
Stretchiness	≥6N
Material	Nitrile
Palm thickness:	min. 0.06mm
Finger thickness:	min. 0.07mm
Sterility	Non-Sterile
Texture	Texture at the fingertips
Internal surface	Powder-Free, Max. 2.0mg/glove
Impermeability level (AQL)	≤1.0
Shape	ambidextrous
Cuff	Rolled
Glove length	≥240mm
Glove palm width	XS: 70mm, S: 80mm, M:95mm, L:110mm, XL: 120mm.
Suitable for sensitive skin	Yes. Latex Free
Shelf Time:	5 years
storage temperature	5-35°C degrees




Product EAN-13 Codes

XS 6 948794 411746	S 6 948794 411753	M 6 948794 411760
L 6 948794 411777	XL 6 948794 411784	XXL 6 948794 411791

Package Specification

Boxes/Carton:	10
Carton G. Weight:	XS-3.9KG S-4.2KG, M-4.5KG, L-4.8KG , XL-5.1KG, XXL-5.8KG
Maximum boxes per EU Pallet:	880 boxes/ pallet
EU Pallets Qty per full truck:	33 pallets
EU Pallet G. Weight:	21.5 kg
G. Weight per pallet for XS	364.7 kg
G.Weight per pallet for S	391.1 kg
G.Weight per pallet for M	417.5 kg
G. Weight per pallet for L	443.9 kg
G. Weight per pallet for XL	470.3 kg
G. Weight per pallet for XXL	531.9 kg



Box size: 240mm x 122mm x 65mm